

**BIO-RAD**  
CFX Opus



Specialist : 蔡鼎諾 (Tim)



Genmall

# qPCR Operation Training



# Outline

01

Principle of PCR

02

Principle of qPCR

03

Precautions of qPCR

04

Operation of qPCR

# PCR

## DNA

## Principles

## Elements

## Procedure

### What is Polymerase Chain Reaction ?

A common laboratory technique which is used to “amplify” - copy - small segments of DNA (In Vitro) via DNA polymerase and a pair of specific primers.





PCR

DNA


Principles

Elements

Procedure



GC pair > AT pair




PCR

DNA

Principles

Elements

Procedure




PCR

DNA

Principles

Elements

Procedure



PCR

DNA

Principles

Elements

Procedure

- ✓ **DNA template** : A segment of DNA which can be used as a model for duplication.
- ✓ **Primer** : Two primers (**forward& reverse primer**) are used in each PCR reaction, and they are designed so that they flank the target region (**region that should be copied**).
- ✓ **DNA polymerase** : Synthesizes - builds - two new strands of DNA, using the original strands as templates.

PCR

DNA

Principles

Elements

Procedure

- ✓ **Deoxynucleotide triphosphate (dNTP)** : There are four types of dNTP 「**dATP**、**dGTP**、**dCTP**、**dTTP**」 in PCR is to expand the growing DNA strand with the help of DNA polymerase. ◦
- ✓ **Magnesium chloride ( $\text{MgCl}_2$ )** : Magnesium ion ( **$\text{Mg}^{2+}$** ) functions as a **cofactor** for activity of DNA polymerases by enabling incorporation of dNTPs during polymerization.
- ✓ **Buffer** : Provides a suitable chemical environment for activity of DNA polymerase. The buffer pH is usually between 8.0 and 9.5 (**KCl**) and is often stabilized by **Tris-HCl**, which promotes primer annealing.


PCR

DNA

Principles

Elements

Procedure





2022 

— q P C R —

...

# qPCR

Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe

1. **Real-Time PCR** : divide into two parts : 「**Real-time**」 and 「**PCR**」
2. Real-Time PCR : It is also known as 「**quantitative PCR**」
3. Abbreviation : 「**Real Time RT-PCR**」 or 「**qPCR**」
4. Method : Add 「**fluorescent DNA binding dye**」 or 「**probe**」 for detection during DNA amplification (**PCR**)
5. Advantage : **qPCR** is sensitive, specific detection (via **fluorescence**) and quantification of nucleic acid targets during PCR cycle (so-called **Real time**) .



# qPCR

Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe

BIO-RAD

## 1 Heat denaturation



## 2 Primer annealing



## 3 Extension



# PCR V.S. qPCR

BIO-RAD

qPCR

Comparison



Principle

$C_T$  value

DNA binding dye

Hydrolysis probe


+Taq  
+dNTP  
+MgCl<sub>2</sub>  
+H<sub>2</sub>O  
+Buffer  
+Template  
+Primer  
(On ice)



- End point
- Fuzzy data
- Labor intensive

3hrs (1.5 hrs labor)

+2x Master Mix  
+H<sub>2</sub>O  
+Template  
+Primer  
(Room Temp)



- Online monitoring
- More sensitive
- Automated

<1 hrs (20 min labor)

qPCR


Comparison

Principle

 $C_T$  value

DNA binding dye

Hydrolysis probe



qPCR

Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe




Fig. 1.1. Amplification plot. Baseline-subtracted fluorescence is shown.

- ✓ **Threshold line** (background or baseline) : At the start of the run, when the amount of PCR product is low, this produces very little fluorescence.
- ✓ **Exponential phase** : Once the amount of PCR product is amplified enough, it will enter the exponential phase. This is when the amount of PCR product doubles for every PCR cycle. [  $C_T$  ] can be used and be defined as quantitative marker during this phase.
- ✓ **Non-exponential Plateau phase** : Once all of the reagents, such as nucleotides, have been used up in the PCR reaction, the amplification will slow and ultimately plateau. This is the region where no more PCR products cannot be produced.



qPCR


Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe



- ✓ **Threshold** : Basic reaction line of **background** which is derived from fluorescent reagents.
- ✓  **$C_T$  (cycle threshold) /  $C_q$  (cycle quantification)** is defined as the number of cycles required for the fluorescent signal to **cross the threshold** (i.e. **exceeds background level**).

qPCR


Comparison

Principle

 $C_T$  value

DNA binding dye

Hydrolysis probe



$C_T$  value are inversely proportional to the **amount** of target nucleic acid

# DNA binding dye – SYBR Green I

BIO-RAD

qPCR

Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe



## 1 Heat denaturation



## 2 Primer annealing



## 3 Extension



SYBR Green I is a dsDNA binding dye, which can be used to quantify amplicon amount during the course of the PCR by tracking overall fluorescence emission.

# DNA binding dye – SYBR Green I

BIO-RAD

qPCR



Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe



How to confirm the 「**specificity**」 of PCR products when performing **SYBR Green** qPCR experiments?

1. Method: Use the characteristics of DNA  **$T_m$  value** to ensure.
2. Analytical method :
  - ① **Melt Curve** ( The upper left )
  - ② **Melt Peak** ( The upper right )



# Melting Curve - T<sub>m</sub> value

BIO-RAD

qPCR


Comparison

Principle

C<sub>T</sub> value

DNA binding dye

Hydrolysis probe



## ➤ Melting temperature (T<sub>m</sub>)

- ① The temperature at which **one half of the DNA duplex** will dissociate to become **single stranded** and indicates the duplex stability. ◦
- ② Depends on nucleotide content ( **GC ratio** ) and **length**

# Melting Curve - T<sub>m</sub> value

BIO-RAD

qPCR


Comparison

Principle

C<sub>T</sub> value

DNA binding dye

Hydrolysis probe



## ➤ Melting temperature (T<sub>m</sub>)

- ① The temperature at which **one half of the DNA duplex** will dissociate to become **single stranded** and indicates the duplex stability. ◦
- ② Depends on nucleotide content (**GC ratio**) and **length**

# Melting Peak – T<sub>m</sub> value

qPCR


Comparison

Principle

C<sub>T</sub> value

DNA binding dye

Hydrolysis probe



Target Product

Primer-dimers

Mispriming

**Fig. 6. Melt curve profiles.** A, a single well-defined peak indicates a single specific product. B, two or more peaks indicate poor specificity. A typical primer-dimer (in red) is due to excess final primer concentrations. Primer-dimers typically amplify in the mid-70°C range. C, two or more peaks indicate poor specificity. A typical mispriming (in red) is shown.

## ➤ Melting temperature (T<sub>m</sub>)

- ① The temperature at which **one half of the DNA duplex** will dissociate to become **single stranded** and indicates the duplex stability. ◦
- ② Depends on nucleotide content (**GC ratio**) and **length**

# Melting Peak – T<sub>m</sub> value

BIO-RAD

qPCR



Comparison

Principle

C<sub>T</sub> value

DNA binding dye

Hydrolysis probe



# Hydrolysis probe – TaqMan probe

BIO-RAD

qPCR


Comparison

Principle

C<sub>T</sub> value

DNA binding dye

Hydrolysis probe



- ① **Structure** : Labeled with a fluorescent **reporter at the 5' end** and a **quencher at the 3' end**.
- ② **Enzyme** : The dsDNA-specific **5' → 3' exonuclease** activity of *Taq*
- ③ **Principle** : The **reporter** is separated from the **quencher**, resulting in a fluorescence signal that is proportional to **the amount of amplified product** in the sample.
- ④ **Advantages** : **High specificity** and the ability to perform **multiplex** reactions.

# Hydrolysis probe – TaqMan probe

BIO-RAD

qPCR

Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe


## Extension Step

1. Strand Displacement

2. Cleavage

3. Polymerization Complete

4. Detection



# Binding dye v.s. Probe

qPCR

Comparison

Principle

$C_T$  value

DNA binding dye

Hydrolysis probe


## ➤ SYBR Green I Dye

- ① Singleplex
- ② Non-specific
- ③ Binding site : **dsDNA**
- ④ Sensitivity : Middle (10 - 100 copies)
- ⑤ Notice : Avoid **Primer-dimer** formation
- ⑥ Cost: **Inexpensive**


## ➤ TaqMan Probe

- ① Multiplex
- ② **Highly Specific**
- ③ Binding site : **probe-specific region**
- ④ Sensitivity : High (1-10 copies)
- ⑤ Notice : **Specific fluorescent probe** is necessary. (Must be less than 30 bp)
- ⑥ Cost: **Expensive**

### A. SYBR Green



### B. TaqMan Probe







**2022** 

Application





# Standard Curve

BIO-RAD

Essential

Absolute Qu.

Relative Qu.

Application

MIQE

Notice

- Linear standard curve : Coefficient of Correlation,  $R^2 > 0.980$
- PCR Efficiency :  $> 90 - 110\%$
- Equation :
  - ① Amplification efficiency,  $E = 10^{-1/\text{slope}}$
  - ②  $\% \text{ Efficiency} = (E - 1) \times 100\%$
- When  $E < 90\%$  : Might be interfered by **PCR inhibitors** in your sample or solution.
- When  $E > 110\%$  : Might be influenced by **pipetting error** or **primer-dimer**.

# Standard Curve

BIO-RAD

Essential


Absolute Qu.

Relative Qu.

Application

MIQE

Notice



Slope: -3.395

Y-Intercept: 40.47

Efficiency: 0.97

R<sup>2</sup>: 0.998

# Essential

## Absolute Qu.

## Relative Qu.

## Application

## MIQE

## Notice

✓ 「Absolute quantification」 is performed by constructing a standard curve for each gene of the genes of interest (GOI).

✓ Plotting the quantification cycle (Cq) values against log[quantity] of a dilution series of known GOI amount. (Interpolation)

✓ According to the different standard products, it is usually expressed in the form of common units as follows:

- ① Concentration unit : ng/ml 、 μg/ml
- ② Measurement unit : ng 、 pg 、 copy number

✓ Note for absolute quantification : The slope of the standard curve and the correlation coefficient ( $R^2$ )

- ① Slope : 「3.0 ~ -3.6」 , (Correspond to PCR efficiency from 90%~110%)
- ② Correlation coefficient : 「 $R^2 > 0.980$ 」

Essential

Absolute Qu.

Relative Qu.

Application


MIQE

Notice

1

$$y = mx + b$$

$$y = -3.50x + 37.20; r^2 = 0.999$$

  $C_T$ 


2

log Quantity

3

Sample	Replicate	$C_T$	Copies
A	1	18.61	204,577
A	2	18.41	234,115
A	3	18.87	172,300
Average			203,664 ± 30,917
B	1	17.06	569,789
B	2	17.07	563,823
B	3	17.00	591,173
Average			574,928 ± 14,381



Essential

Absolute Qu.

Relative Qu.

Application

MIQE

Notice

	Reference	Target
Tissue #1 (Control):	21	22
Tissue #2 (Test) :	20	24
<hr/>		
1 <sup>st</sup> Delta	$\Delta C_T \#1:$	$22 - 21 = 1$
	$\Delta C_T \#2:$	$24 - 20 = 4$
<hr/>		
2 <sup>nd</sup> Delta	$\Delta\Delta C_T :$	$1 - 4 = -3$

$$\text{Fold induction} = 2^{-\Delta\Delta C_T} = 2^{-(-3)} = 8$$

Essential

Absolute Qu.

Relative Qu.

Application

MIQE

Notice

### Gene Expression

- ① Drug reaction
- ② Tumor marker
- ③ Gene regulation (siRNA/miRNA)
- ④ Gene therapy
- ⑤ Microarray
- ⑥ GMO detection

### Pathogen detection

- ① Multiplex detection
- ② Multiple specimen detection
- ③ Synchronous analysis : Qualitative / Quantitative
- ④ Drug Safety and Effectiveness Monitoring

### SNP Genotyping Analysis

- ① Evaluate the identified position of SNP
- ② Use for drug design
- ③ Genetic disease detection

Essential

Absolute Qu.

Relative Qu.

Application

MIQE

Notice




BIO-RAD

Clinical Chemistry 55:4  
611–622 (2009)

Special Report

## The MIQE Guidelines: Minimum Information for Publication of Quantitative Real-Time PCR Experiments

Stephen A. Bustin,<sup>1\*</sup> Vladimir Benes,<sup>2</sup> Jeremy A. Garson,<sup>3,4</sup> Jan Hellemans,<sup>5</sup> Jim Huggett,<sup>6</sup>  
Mikael Kubista,<sup>7,8</sup> Reinhold Mueller,<sup>9</sup> Tania Nolan,<sup>10</sup> Michael W. Pfaffl,<sup>11</sup> Gregory L. Shipley,<sup>12</sup>  
Jo Vandesompele,<sup>5</sup> and Carl T. Wittwer<sup>13,14</sup>





# Essential


# Absolute Qu.

# Relative Qu.

# Application

# MIQE

# Notice







# MIQE guideline

1

檢核項目	重要性	說明
<b>EXPERIMENTAL DESIGN</b>		
Definition of experimental and control groups	E	敘述實驗設計組別、樣本數以及執行的實驗室等
Number within each group	E	
Assay carried out by core lab or investigator's lab?	D	
Acknowledgement of authors' contributions	D	

1

<b>SAMPLE</b>		
Description	E	敘述樣本來源(例如組織或細胞)、細胞數(量)、處理方式(多少時間內完成)、保存環境(溫度)
Volume/mass of sample processed	D	
Microdissection or macrodissection	E	
Processing procedure	E	
If frozen - how and how quickly?	E	
If fixed - with what, how quickly?	E	
Sample storage conditions and duration (especially for FFPE samples)	E	

2 3

<b>NUCLEIC ACID EXTRACTION</b>		
Procedure and/or instrumentation	E	提供萃取的方法或儀器
Name of kit and details of any modifications	E	提供 <b>Kit</b> 的名稱/型號/批號/廠牌
Source of additional reagents used	D	是否有額外添加的試劑，有則註明型號廠牌
Details of DNase or RNase treatment	E	敘述對 <b>DNase (for RNA 實驗)</b> or <b>RNase (for DNA 實驗)</b> 的處理
Contamination assessment (DNA or RNA)	E	以尚未進行反轉錄的 <b>RNA</b> 樣本進行 <b>qPCR</b> ，確認是否有產物
Nucleic acid quantification	E	核酸定量:使用的儀器廠牌/方法， <b>OD260/OD280</b> 比值，
Instrument and method	E	計算濃度和總量
Purity (A260/A280)	D	確認 <b>RNA</b> 品質: 求 <b>RIN/RQI</b> 值 <sup>1</sup> ，跑膠確認 <b>RNA</b> 是否降解以及測試樣品內是否有抑制物質(效應) <sup>2</sup>
Yield	D	
RNA integrity method/instrument	E	
RIN/RQI or Cq of 3' and 5' transcripts	E	
Electrophoresis traces	D	
Inhibition testing (Cq dilutions, spike or other)	E	

4

<b>REVERSE TRANSCRIPTION</b>		
Complete reaction conditions	E	敘述反轉錄反應條件
Amount of RNA and reaction volume	E	<b>RNA</b> 量與反應體積
Priming oligonucleotide (if using GSP) and concentration	E	<b>Primer</b> 種類(廠牌)和濃度
Reverse transcriptase and concentration	E	反轉錄酵素(廠牌)和濃度
Temperature and time	E	溫度/時間/cycle 數
Manufacturer of reagents and catalogue numbers	D	試劑廠牌和型號
Cqs with and without RT	D	<b>RT</b> 試劑測試
Storage conditions of cDNA	D	記錄保存條件



# MIQE guideline

5

qPCR TARGET INFORMATION		
If multiplex, efficiency and LOD of each assay.	E	於同一試管做多基因定量才需要確定該項目
Sequence accession number	E	先以資料庫或軟體確定基因序列、比對、二級結構分析、Primer 對應於 exon 位置等資訊
Location of amplicon	D	
Amplicon length	E	
<i>In silico</i> specificity screen (BLAST, etc)	E	
Pseudogenes, retropseudogenes or other homologs?	D	
Sequence alignment	D	
Secondary structure analysis of amplicon	D	
Location of each primer by exon or intron (if applicable)	E	
What splice variants are targeted?	E	

5

qPCR OLIGONUCLEOTIDES		
Primer sequences	E	Primer 或 Probe 合成的序列、廠商、純化方法等
RTPrimerDB Identification Number	D	
Probe sequences	D <sup>3</sup>	
Location and identity of any modifications	E	
Manufacturer of oligonucleotides	D	
Purification method	D	

5

qPCR PROTOCOL		
Complete reaction conditions	E	敘述定量 PCR 實驗反應條件
Reaction volume and amount of cDNA/DNA	E	cDNA 量與反應體積
Primer, (probe), Mg <sup>++</sup> and dNTP concentrations	E	Primer、鎂離子濃度、聚合酶濃度、緩衝液組成、其他添加物等
Polymerase identity and concentration	E	
Buffer/kit identity and manufacturer	E	
Exact chemical constitution of the buffer	D	
Additives (SYBR Green I, DMSO, etc.)	E	分析盤(管)廠牌和型號
Manufacturer of plates/tubes and catalog number	D	
Complete thermocycling parameters	E	反應條件(溫度/時間/cycles)
Reaction setup (manual/robotic)	D	反應設定
Manufacturer of qPCR instrument	E	儀器廠牌機型

6

qPCR VALIDATION		
Evidence of optimisation (from gradients)	D	是否有測試 annealing 溫度梯度
Specificity (gel, sequence, melt, or digest)	E	確認反應專一性
For SYBR Green I, Cq of the NTC	E	是否 NTC 控制組有產物
Standard curves with slope and y-intercept	E	建立標準曲線 <sup>4</sup> :由斜率推算 PCR 效率、計算信賴區間與標準差以及 r <sup>2</sup>
PCR efficiency calculated from slope	E	
Confidence interval for PCR efficiency or standard error	D	
r <sup>2</sup> of standard curve	E	
Linear dynamic range	E	求得線性關係範圍 <sup>5</sup> : 計算信賴區間、最低極限時的 Cq 變異數以及是否確認 LOD 偵測極限
Cq variation at lower limit	E	
Confidence intervals throughout range	D	
Evidence for limit of detection	E	於同一試管做多基因定量計算效率和偵測極限
If multiplex, efficiency and LOD of each assay.	E	

DATA ANALYSIS		
qPCR analysis program (source, version)	E	分析程式/工具(版本)、分析方法 (threshold 或 regression)以及敘述樣本盤排列方式等
Cq method determination	E	NTC 控制組是否有產物
Outlier identification and disposition	E	
Results of NTCs	E	Reference gene 的數目和選擇
Justification of number and choice of reference genes	E	敘述標準化方法
Description of normalisation method	E	生物樣本數
Number and concordance of biological replicates	D	單一生物樣本在儀器偵測的 n 值
Number and stage (RT or qPCR) of technical replicates	E	
Repeatability (intra-assay variation)	E	計算重複性、再現性、使用何種統計方法、統計軟體(版本)
Reproducibility (inter-assay variation, %CV)	D	
Power analysis	D	
Statistical methods for result significance	E	
Software (source, version)	E	以 RDML 格式提供原始數據
Cq or raw data submission using RDML	D	

Essential



Absolute Qu.

Relative Qu.

Application

MIQE

Notice

**Biological Replicates****Technical Replicates**

**Fig. 6. Experimental replicates.** All experiments should be designed with a combination of biological and technical replicates. This illustrates a simple experiment with triplicate biological samples from control (■) and treatment/experimental (■) conditions. For each biological sample, three technical replicates are recommended for the gene of interest as well as for the reference gene(s). This results in a total of at least 36 samples plus the duplicate NTC (■) for a total of >40 wells.

# Common 「PCR inhibitor」

BIO-RAD

Essential

Absolute Qu.

Relative Qu.

Application

MIQE

Notice



Sample	Assay
Bile salts	Cholic acid Deoxycholic acid
Polysaccharides	Polygalacturonic acid
Lipids	Cholesterol hydrochloride
Algae	Alginic acid
Clay	Montmorillonite
Blood	Heparin Hematin EDTA Serum
Soil	Humic acid
Textile	Indigo
Wine	Tannic acid
Plants	Cellulose Pectin (for fiber control)
Hair, tissues	Melanin
Bones, teeth	CaCl <sub>2</sub>

Sample	Assay
Cell lysates	MEM+FBS PBS Trypsin
Sample preparation	NaAc NaCl EtOH Isopropanol TRIzol InstaGene Matrix
Miscellaneous	Green tea Chocolate Dust
	SDS DMSO DTT
Body fluids	Spermidine Urea

★Red : From 「sample」

★Green : From 「reagent / solution」



# Amplicon

Essential

Absolute Qu.


Relative Qu.

Application

MIQE

Notice

1. Recommended length : 75 - 200bp
2. Regional structure : Should be free of strong secondary structures .
3. Repeats : Prevent quadruple repeats of G or C appear in designed region.
4. CG % : 40 - 60%



# Primer design

Essential

Absolute Qu.


Relative Qu.

Application

MIQE

Notice

1. Recommended length : 18 - 25bp
2. Regional structure : Should be free of strong secondary structures and self-complementarity 「Primer-dimer」
3. Repeats : Prevent triple repeats of G or C appear in designed region.
4. CG % : 40 - 60%
5. T<sub>m</sub> value : 50 - 65°C , Formula :  $T_m = 4^{\circ}\text{C} (G + C) + 2^{\circ}\text{C} (A + T)$



Taq DNA polymerase

DNA

# (q)PCR - 產線支援

BIO-RAD

Essential

Absolute Qu.

Relative Qu.

Application


MIQE

Notice

PCR



CFX Manager Software



(q)PCR Reagents



Real-Time PCR (qPCR)



(q)PCR consumables



Automation options





**2022** 

— Operation —





# System Overview - Front View

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# System Overview - Back View

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Touch Screen Overview - Open lid

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Creating a Protocol

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Modifying the Settings in a Protocol Step

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Modifying the Settings in a Protocol Step

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Modifying the Settings in a Protocol Step

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics





# Changing Target Gradient and Ramp Rate

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Changing Target Gradient and Ramp Rate

BIO-RAD

Appearance

Device Boot

Consumables

Software

Optics

### Step Options

Temperature:  °C

Time:  HH:MM:SS

Ramp Rate:  °C/s

Increment:  °C/cycle

Extend:  s/cycle

Beep: ☐

Plate Read: ☐

OK

# Adding or Removing a Temperature Gradient

BIO-RAD

Appearance

Device Boot

Consumables

Software


Optics


### Step Options


Temperature:

Gradient (°C): ☒

Time:  HH:MM:SS

Ramp Rate:  °C/s 

Increment:  °C/cycle 

Extend:  s/cycle 



Beep: ☐ 

Plate Read: ☐ 

OK

A	100
B	99.8
C	99.2
D	98.2
E	97.1
F	96.2
G	95.5
H	95

# Saving a Protocol

BIO-RAD

Appearance

Device Boot

Consumables

Software

Optics

New Protocol

CarlN

Name:  Volume:   $\mu\text{l}$  Lid Temp:   $^{\circ}\text{C}$

1	2	3	4
95 $^{\circ}\text{C}$	95 $^{\circ}\text{C}$	55 $^{\circ}\text{C}$	
3:00	0:10	0:30	GOTO Step 2 39 X

Save Next

# Saving a Protocol



Appearance

Device Boot

Consumables

Software

Optics



# Run Setup

BIO-RAD

Appearance

Device Boot

Consumables

Software

Optics

New Protocol

CarlN

Name:  Volume:   $\mu\text{l}$  Lid Temp:   $^{\circ}\text{C}$

1	2	3	4
95 $^{\circ}\text{C}$	95 $^{\circ}\text{C}$	55 $^{\circ}\text{C}$	
3:00	0:10	0:30	GOTO Step 2 39 X

Insert Remove Options Save **Next**

Appearance

Device Boot

Consumables

Software

Optics

# Run Setup - Start Run

BIO-RAD


[Back](#) Run Setup CARLN



Name: 3AMPMelt\_CarLN Volume:   $\mu$ l Lid Temp:   $^{\circ}$ C

Scan Mode: ☒ SYBR/FAM ☐ All Channels ☐ FRET

Plate ID:

Run File Name:

Save Location:  CARLN\...\CarLN

Notification: ☒  ☒  cnavar@celltech.com



# Touch Screen Overview - Proceeding

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Touch Screen Overview - File

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# The File Browser Screen

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Extracting and Editing a Protocol from a Run

BIO-RAD

Appearance

Device Boot

Consumables

Software

Optics

Run Setup

Name: 3AMPMelt\_CarlN Volume: 30  $\mu$ l Lid Temp: 105  $^{\circ}$ C

Scan Mode: ☒ SYBR/FAM ☐ All Channels ☐ FRET


Plate ID:

Run File Name: 3AMPMelt\_CarlN\_20191117\_131432\_OPUS0001\_CARLN

Save Location: CARLN\...\CarlN

Notification: ☒ ☒ cnavar@celltech.com

Open Lid Run



# Consumables - Specification

BIO-RAD

Appearance

Device Boot



Consumables

Software


Optics



Low – Profile Only !!



Well Color Effect on Signal Strength





# Notice for usage

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Connection with Opus 96 - Overview

BIO-RAD


Appearance

Device Boot


Consumables

Software

Optics






BIO RAD CFX Maestro



Protocol

Plate

Start Run



Start Run !



# Operation setting - Overview

BIO-RAD


Appearance

Device Boot

Consumables

Software

Optics



# Optical design for detection

BIO-RAD




Appearance

Device Boot

Consumables

Software

Optics



# Optical design for detection

BIO-RAD


Appearance

Device Boot


Consumables

Software


Optics



BIO-RAD



Competitor



Vertical light path : 「Rox dye」 calibration is unnecessary

Fluorescence : 「high strength / specificity」

Light source : **LEDs** with long lifespan, allowing you to use them right away

# Optical application - Dye / Wave band

BIO-RAD



Appearance

Device Boot

Consumables

Software

Optics



Channel	Excitation (nm)	Detection (nm)	Calibrated Fluorophores
1	450-490	515-530	FAM™, SYBR Green I™, EvaGreen™
2	515-535	560-580	VIC®, HEX™, TET™, Cal Gold 540™
3	560-590	610-650	ROX™, Texas Red®, Cal Red 610™
4	620-650	675-690	Cy5, Quasar 670™
5	672-684	705-730	Quasar 705™
6	450-490	560-580	Accommodates <b>FRET Chemistry</b>



**2022** 

— Supplement —







Blocks

Gradient

Application



Mass-reduced sample block



Uniform / Rapid heating



Max ramp rate	5°C / sec
Temp Accuracy	± 0.2°C
Temp Uniformity	± 0.3°C

Accurate temperature control

BIO-RAD



Uniformity ↑




Supreme 「 uniformity 」 in any position


Blocks

Gradient

Application




Graphical interface : friendly operation and setting



Peltier Thermal System

		1	2	3	4	5	6	7	8	9	10	11	12
A	70.0	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
B	69.5	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
C	68.4	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
D	66.4	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
E	64.0	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
F	62.0	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
G	60.7	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk
H	60.0	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk	Unk

Gradient temperature setting : 「1 - 24°C」



Same starting point / Reaction is accurate




30s



# Blocks

# Gradient

# Application



**Optimization of an assay results in better yields and specificity.** Results show that assays A and B can be run at an **annealing temperature** of **57.9°C** on the same plate. Higher temperatures result in a reduced yield in assay A while lower temperatures result in nonspecific products in assay B. M, markers.



**正茂生物科技股份有限公司**  
Genmall Biotechnology Co., Ltd.



*Thank You!*

