


CPR + AED

心肺復甦術程序




AED 簡單操作四步驟 AED STEPS

- 1 打開電源**
Turn on the power.
- 2 貼上電擊片**
Attach pads.
- 3 分析心律**
此時不可接觸或移動病人
Analyze heart rhythm
(make sure no one touches the victim)
- 4 電擊**
若機器有指示需要電擊
Shock (if AED advises to shock the person)



Attention:

- 1** 電擊貼片黏貼位置如有潮溼，必須擦乾，並緊貼於皮膚(如有外傷亦需避開)。
Wipe the chest dry if it is wet, attach pads to bare chest. (Avoid attaching on wound)
- 2** 緊急醫療救護法第14-2條：救護人員以外之人使用AED救人免責(救護人員於非值勤期間，亦適用之)。
Emergency Medical Services Act: People other than paramedics use AED do not have any responsibilities.



電擊貼片貼附注意事項

- 金屬配飾需先移除
Remove all metal jewelries.
- 患者上半身及地面須保持乾燥
Keep body and the ground dry.
- 電擊貼片需緊貼皮膚，避開傷口或衣物、貼布
Attach pads to patient.