

Applicable reagents

Luminescence

ECL™	ECL Plus™	ECL Advance™
Lumi-Light Plus	SuperSignal™	CDP-Star™
CSPD™	Immobilon	Bright-Star™

Fluorescence/IR

DY676	DY781
Alexa Fluor™ 680	Alexa Fluor™ 700

Fluorescence/red

Alexa Fluor™ 633	Alexa Fluor™ 635	Alexa Fluor™ 647
Cy™ 5	BODIPY™ 650/665	DiD
TOTO™ -3	DDAO phosphate	

Fluorescence/green

SYPRO™ Red	Cy™ 3	TAMRA™
ROX™	HEX™	Alexa Fluor™ 532
Alexa Fluor™ 546	Deep Purple	Pro-Q™ Diamond
Rhodamine Red™	BODIPY™ 576/589	NED™
R-phycoerythrin	RFP	HNPP
Tetramethylrhodamine		

Fluorescence/blue

SYBR™ Green I	SYBR™ Green II	SYBR™ Gold
SYPRO™ Ruby	SYPRO™ Orange	SYPRO™ Tangerine
FITC	FAM™	EGFP
ECFP	AttoPhos™	

Fluorescence/epi-UV

Ethidium Bromide	SYPRO™ Rose	Qdot™
------------------	-------------	-------

Fluorescence/trans-UV

Ethidium Bromide

Digitizing/white epi-illumination/white transillumination

Silver staining	CBB	NBT/BCIP
X-ray film, etc.		

* The reagents listed above are representative reagents with high performance results. Please contact us for other reagents.



FUJIFILM 台灣總代理／
恒昶實業公司
台北市民權東路六段38號 電話：(02) 2791-1188 (總機)